

Tackle Regulatory Challenges

NanoPhotometer® CFR21 Software

Compliance for Electronic Records



Electronic Signature

Save measurement data by confirming login name and password. Saved files provide user name/author, User ID, date and time.



User Management

Facilitate access of shared data and stored methods by creating multiple user accounts with different access rights. User role options are Administrator, Power User, and User.



Audit Trail

All actions and preference changes are recorded in an audit log which contains a sequence ID, time stamp, User ID, and category for each action.



Data Integrity

Measurement data are saved as a proprietary Hash Code secured and encrypted IDS file and can also be stored in PDF format.



Data Backup

A proprietary IDS file is automatically saved as a backup file in an Autosave folder accessible by Administrators. Entire system backups can be performed and saved as a copy on a USB flash drive.





Access Control

Device access is limited to authorized individuals with password protected user accounts. Accessibility to folders and data depends on user's access rights.

The CFR21 software complies with FDA 21 CFR Part 11 requirements and is an optional software tool ideal for GxP laboratories, which require proper electronic record keeping. It includes user management, access control, electronic signatures, data integrity, security, and audit trail functionality.

User Management

Individual Role Based Access Control (RBAC) provides password protected access and control of the NanoPhotometer®. Create multiple user accounts with different access rights which are handled in a hierarchic structure. User role options are Administrator, Power User, and User. Organize users into working groups to facilitate access of shared data and stored methods within a lab. Option for increased transparency with Four Eye Administrator Rights. Various password settings are available within the CFR21 software – for example secure password and password expiration options. Effectively improve data security and fulfill audit regulations easily with flexible and appropriate RBAC user management solutions. All features can be enabled or disabled on demand to meet your laboratory needs.



Electronic Signature

Measurement data are saved by confirming login name and password. Saved files include the user name/author, login name, date and time for proper electronic record keeping. IDS and PDF files cannot be altered and ensure data integrity.

Implen NanoPhotometer®

Instrument Type	NP80		
Version	NPOS 4.2 build 14756		
Serial Number	M80798		
Selftest passed	2019-09-05; 13:02		
Autosave	No		
File Name	GroupA/rmsmithBSA.lds		
Reason	Read/Save/Print		
User ID	bjones		
User Name	Mark Smith	Becky Jones	
eSign Date	2019-09-05	2019-09-05	
eSign Time	13:29:46	13:39:16	

Parameter	Protein UV	Wavelength (nm)	280
Method	BSA	Background Correction	320 nm
Type	NanoVolume	Air Bubble Recognition	Off
Mode	1.500	Manual Dilution Factor	1.000
Protein Factor	1-2		
Volume (ul)			

#	Name	Conc.	Units	A230	A260	A280	A320	A260/A280	Dilution
1	Blank 1	0.0000	mg/ml	0.000	0.000	0.000	0.000	0.000	
2	Sample 1	0.0000	mg/ml	0.002	0.001	0.006	0.006	---	15
3	Sample 2	3.3345	mg/ml	18.02	2.756	2.230	0.007	1.237	15
4	Sample 3	3.3435	mg/ml	18.04	2.752	2.234	0.005	1.232	15

Audit Trail

The audit trail automatically records all actions and preference changes in an audit log. The audit log contains a sequence ID, time stamp, user ID and category for each action. Audit trails are sorted by year and can be navigated with a search function. For documentation purpose the audit trail can be saved in PDF format by an Administrator.

Audit Trail

ID	Date/Time	Login Name	Category	Action	Details
05462	2022-03-15 17:38:27	Admin	Administrator	Sample measurement	109147.1
05461	2022-03-15 17:38:18	Admin	Administrator	Sample measurement	109147.2
05460	2022-03-15 17:38:12	Admin	Administrator	Sample measurement	109147.1
05459	2022-03-15 17:38:06	Admin	Administrator	Sample measurement	109147
05458	2022-03-15 17:37:54	Admin	Administrator	Sample measurement	109150. Item 2
05457	2022-03-15 17:37:39	Admin	Administrator	Sample measurement	109150. Item 1
05456	2022-03-15 17:37:26	Admin	Administrator	Sample measurement	109150. Item
05455	2022-03-15 17:37:07	Admin	Administrator	Sample measurement	109150.2
05454	2022-03-15 17:36:58	Admin	Administrator	Sample measurement	109150.1
05453	2022-03-15 17:36:52	Admin	Administrator	Sample measurement	109150
05452	2022-03-15 17:36:45	Admin	Administrator	Sample measurement	Air
05451	2022-03-15 17:36:45	Admin	Administrator	Blank measurement	Blank 2
05450	2022-03-15 17:36:02	Admin	Administrator	Blank measurement	Blank 1
05449	2022-03-15 17:31:43	Admin	Administrator	Login	NanoPhotometer
05448	2022-03-15 17:04:36	System	System	Warning message	Battery is low. Please connect to power!
05447	2022-03-15 17:13:50	Admin	Administrator	Logout	Control Device
05446	2022-03-15 17:10:31	Admin	Administrator	File opened successfully	NanoPhotometer/Admin/Control_3M.lds
05445	2022-03-15 12:10:17	Admin	Administrator	Login	Control Device

Features and specifications are subject to change without notice.