


## Validating Absorbance Accuracy and Reproducibility with Neutral Density Filters on the Implen OD600® Device



For Validation using the Neutral Density Glass filters, open the OD600 app on the Implen OD600® device.

### Step 1: Validate Absorbance Accuracy at 1A with 666-F4

Blank against air.

Put the first certified Neutral Density Glass filter (**666-F4**) into the cuvette port. Light path is back to front. Place filter in correct position.

Press sample.

Compare results to the specified absorbance in the calibration certificate.

Tolerance is  $\pm 0.01$  A.

### Step 2: Validate Absorbance Accuracy at 3A with 666-F303

Blank against air.

Put the second certified Neutral Density Glass filter (**666-F303**) into the cuvette port. Light path is back to front. Place filter in correct position.

Press sample.

Compare results to the specified absorbance in the calibration certificate.

Tolerance is  $\pm 0.03$  A.

Standard	Ident-No.	Ordinate Reading (Absorbance) $\pm$ MU(*) at the following wavelengths:				
		440 nm	465 nm	546 nm	590 nm	635 nm
666-F4	102950	1.010 $\pm 0.010$	0.938 $\pm 0.010$	0.950 $\pm 0.010$	1.016 $\pm 0.010$	0.984 $\pm 0.010$
666-F303	E0102	3.305 $\pm 0.030$	3.063 $\pm 0.030$	3.005 $\pm 0.030$	3.076 $\pm 0.030$	2.891 $\pm 0.030$

(\*) MU: Measurement Uncertainty

### Step 3: Validate Reproducibility

Blank against air.

Repeat 4 measurements with the 666-F4 filter. Tolerance is  $\pm 0.01$  A.

### Step 4: Validate Baseline

Blank against air.

Repeat 4 measurements with air.

Tolerance is  $\pm 0.002$  A.



*If there are any problems or questions with the IQ/OQ documentation please do not hesitate to contact the Implen Support Team: [support@implen.de](mailto:support@implen.de) / +49-89-7263718-20.*