Implen **OD600**[®]

User Manual

Version 1.6

Software Version 1.0.5

The end user of the Implen OD600® product ("End User") hereby takes full responsibility for safe storage and backup of all files and/or data that may be created, saved on or transferred from the device. End User acknowledges that it is possible that data and/or files may be lost or damaged, and further acknowledges and agrees that it has sole responsibility to maintain all appropriate backup of files and data. By using the Implen OD600® device, End User hereby agrees to these terms, and agrees that Implen shall not be held liable for any loss, deletion or damage of any data or files for any reason, including any damages attributable thereto.

Telephone support is available using one of the following phone numbers from your geographic region:

**Europe, Asia, South Pacific,
Middle East, Africa**

Phone: +49 89 72637180
Fax: +49 89 726371854
Email: info@implen.de
Website: www.implen.de

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Implen, Inc.
Unit 104
31194 La Baya Drive
Westlake Village, CA 91362
USA

Declaration of Conformity

1. *Product model / product:*
Product Implen OD600®
2. *Manufacturer / Authorized Representative:*
Name Implen GmbH
Address Schatzbogen 52, D-81829 Munich, Germany
3. *This declaration of conformity is issued under the sole responsibility of the manufacturer.*
4. *Object of the declaration:*
Product Implen OD600®
5. *The object of the declaration described above is in conformity with the relevant Union harmonization legislation:*
2014/35/EU Low Voltage Directive
2014/53/EU Radio Equipment Directive
2014/30/EU The Electromagnetic Compatibility Directive
2011/65/EU The Restriction of Hazardous Substances Directive
FCC 47 CFR Part15 §15.107 and §15.109
6. *References to the relevant harmonized standards used or references to the other technical specifications in relation to which conformity is declared:*

Reference & Date	Title
ETSI EN 301 489-1 V2.2.3	Radio and ancillary equipment for portable use (portable equipment); EUT Operating frequency range: 2.4 - 2.4835 GHz
ETSI EN 300 328 V2.2.2	Electromagnetic compatibility and Radio spectrum Matters (ERM); Wideband transmission systems; Data transmission equipment operating in the 2,4 GHz ISM band and using wide band modulation techniques; Harmonized EN covering the essential requirements of article 3.2 of the R&TTE Directive
ETSI EN 301 489-17 V3.2.4	Electromagnetic compatibility and Radio Spectrum Matters (ERM)
UN38.3, 2013/56/EU	Battery certification and transport test
IEC 61010-1:2010	Safety requirements for electrical equipment for measurement, control, and laboratory use. General requirements.
EN IEC 61326-1:2021	Electromagnetic compatibility- generic emission standard electrical equipment for measurement, control and laboratory use.
EN 18031-1:2024	Common security requirements for radio equipment
IEC 60529:1989/AMD2:2013	Degrees of protection prohibited by enclosures (IP 20)

The product is SGS certified in accordance with UL 61010-1, 3rd Edition (Rev. July 19, 2019) and CAN/CSA C22.2 No. 61010-1, 3rd Edition (A1:2018).

7. *Additional information:*
The technical documentation for the machinery is available from:

Name: Implen GmbH
Address: Schatzbogen 52, D-81829 Munich, Germany
Signed for and behalf of: Implen GmbH
Place of issue: Munich, Germany
Date of issue: 01.08.2025
Name: Martin Sahiri
Function: Managing Director
Signature: 

Contents

1. Implen OD600® at a glance	8
Overview	8
Compatible sample container types	8
Basic Layout	9
Instrument Rear Panel	9
Instrument Bottom View	10
Accessories	10
Standard Accessories	10
Optional Accessories	10
Connectivity	11
Implen OD600® Specifications	12
2. Getting Started	13
Implen OD600® Installation	13
3. Implen OD600® Basics	15
Applications Overview	15
General Features	17
Sample List	17
Sample Handling Tips	18
Operating Conditions	19
Power Management Features	19
Battery Operation	19
Auto-Shutdown Features	19
4. Implen OD600® Applications	20
OD600 Application	20
Method Overview	20
Measurement Protocol	21
Calculations	22
McFarland	23
Method Overview	23
Container selection & Calibration	23
Measurement Protocol	27
Calculations	29
DiluCell™ Technology	30
Method Overview	30
Built-in Virtual Dilution	30
Operating Instructions	30
Specifications	31
5. Preferences	32
Main Menu Options	32
Settings	32
Storage	41
About	41
Help	41
6. Troubleshooting	42
Self-test	42

7. Important Error Messages	43
8. Important Warning messages	44
System Freeze	44
9. Maintenance	45
Maintenance-free technology	45
Software update	45
Contact Implen	45
Replacement Parts	46
Cleaning and General Care	46
De-contamination Procedure	46
10. Warranty	46
11. Appendix 1 - Instalace Implen OD600®	47
Bezpečnostní informace	47
Vybalení a umístění	47
12. Appendix 2 - Implen OD600® Installatie	48
Veiligheidsinformatie	48
Uitpakken en plaatsen	48
13. Appendix 3 - Implen OD600® Installation	49
Informations de sécurité	49
Déballage et mise en place	50
14. Appendix 4 - Implen OD600® Installation	51
Sicherheitsinformationen	51
Auspacken und Platzieren	51
15. Appendix 5 - Installazione di Implen OD600®	52
Informazioni di sicurezza	52
Disimballaggio e collocazione	52
16. Appendix 6 - Instalacja Implen OD600®	53
Informacje dotyczące bezpieczeństwa	53
Rozpakowywanie i umiejscowienie	54
17. Appendix 7 - Implen OD600® Instalación	55
Información de seguridad	55
Desembalaje y colocación	55
18. Appendix 8 - Implen OD600® Installation	56
Säkerhetsinformation	56
Uppackning och placering	56
19. Appendix 9 - Instalação do Implen OD600®	57
Informações de Segurança	57
Desembalagem e posicionamento	58

20. Appendix 10- Инсталиране на Implen OD600®	59
Информация за безопасност.....	59
Разопаковане и поставяне	59
21. Appendix 11 – Implen OD600® Installation	60
Sikkerhedsinformation	60
Udpakning og placering	60
22. Appendix 12 – Implen OD600® Paigaldus	61
Ohutusala teave.....	61
Pakendi lahtipakkimine ja seadme asukoha määramine	62
23. Appendix 13 – Implen OD600® Asennus	62
Turvaohjeet.....	62
Purkaminen ja sijoittaminen	63
24. Appendix 14 – Εγκατάσταση Implen OD600®.....	64
Πληροφορίες ασφάλειας	64
Αποσυσκευασία και τοποθέτηση.....	64
25. Appendix 15 – Implen OD600® Telepítés.....	65
Biztonsági információk	65
Kicsomagolás és elhelyezés	66
26. Appendix 16 – Implen OD600® Instalēšana	67
Drošības informācija	67
Iepakojuma atvēršana un novietošana	67
27. Appendix 17 – Implen OD600® Diegimas	68
Saugaus Naudojimo Informacija	68
Išpakavimas ir Įrengimas	69
28. Appendix 18 – Implen OD600® Installazzjoni.....	69
Informazzjoni dwar is-Sigurtà.....	69
Ippakkjar u Pożizzjoni	70
29. Appendix 19 – Implen OD600® Instalare	71
Informații de Securitate	71
Despachetare și Poziționare	71
30. Appendix 20 – Implen OD600® Inštalácia.....	72
Bezpečnostné Informácie.....	72
Rozbaľovanie a Umiestnenie	73
31. Appendix 21 – Implen OD600® Namestitev	73
Varnostne Informacije.....	73
Razpakiranje in Postavitev.....	74
32. Appendix 22 – Implen OD600® Instalacija.....	75
Sigurnosne informacije	75
Raspakiranje i postavljanje	75

33. Appendix 23 – Implen Suiteáil Implen OD600®.....	76
Eolas Sábháilteachta	76
Díphacáil agus Suíomh.....	76
34. Appendix 24 – Implen OD600®安装.....	77
安全信息.....	77
•开箱和配置.....	77
35. Appendix 25 – Implen OD600® のインストール.....	78
安全情報.....	78
開梱と配置.....	79

1. Implen OD600® at a glance

Overview

Introducing your Implen OD600®, a small and mobile device designed for ease of use in various laboratory settings. Tailored for measuring growth rates of all types of cell, such as E. coli and yeast, this compact instrument offers unparalleled versatility and precision.

Featuring a glove-compatible touchscreen display, it offers direct access to pre-programmed OD600 and McFarland applications.

The Implen OD600® accommodates a range of sample containers without the need for adapters, including DiluCells™, cuvettes, and test/falcon tubes.

With an integrated battery lasting up to 8 hours, uninterrupted operation is ensured. Data export is simple via USB in .csv and .xls formats, and advanced connectivity options like LAN/Ethernet and Wi-Fi® for seamless LIMS integration and automatic time adaptation.

The system offers lifetime accuracy and precision without the need for maintenance or recalibration.

Compatible sample container types

Please note that the light beam is directed from the front of the instrument to the back (as indicated by the arrows on the device); therefore please ensure the cell is inserted in the correct alignment.

Note: The optical height of the device is 8.5 mm.

Note: Tube Bottom Shape: For accurate readings, it is essential that the tube is stable at the optical height of 8.5 mm. Tubes with a round or flat bottom allow the tube to sit securely at this height, ensuring a stable and correct reading. Conical or pointed-bottom tubes are not suitable because they do not align well with the optical path.

Note: Tube transparency: sample containers must be fully transparent.


Table 1: The sample compartment is compatible with the following container types.

Container Type	Path length	Min. Volume	Virtual Dilution	Adapter
DiluCell™ 10 (8.5 mm center height)	1 mm	300 µl	1:10	No
Macro cuvette (max fill volume 4.5 ml)	10 mm	1 ml	No	No
Semi-micro cuvette (max fill volume 1.4 ml)	10 mm	0.5 ml	No	No
10 mm diameter tubes	n/a*	0.9 ml	No	No
12 mm diameter tubes	n/a*	1.1 ml	No	No
16 mm diameter tubes	n/a*	1.8 ml	No	No
18 mm diameter tubes	n/a*	2.0 ml	No	No


*Depends on tube being used.

IMPORTANT WARNING: Always wear protective clothing when handling bacteria or other cells.

Basic Layout

Instrument Rear Panel

Instrument Bottom View

Model name, device serial number and FCC ID are located on the identification plate on the bottom of the instrument.

Accessories

Standard Accessories

- **Charging Cable with Adapter**

The USB-C cable provided makes it possible to charge the Implen OD600® in a normal power outlet. Use the power adapter to plug in the device for use and/or to charge the battery.

Note: Use only the power adapter, supplied with your instrument or a replacement part from the manufacturer or your supplier.

Optional Accessories

- **Neutral Density Glass Filter**

The Implen OD600® is recalibration free; therefore, it is not necessary to control the photometric accuracy on a regular basis. If the SOPs of a laboratory require a routine control of the photometric accuracy, a certified neutral density glass filter can be used as a control of photometric accuracy for your OD600® applications.

- **DiluCell™**

Description

The DiluCell™ is a disposable cuvette with a shortened path length for virtual dilution of cuvette-based samples.

The DiluCell™ allows an automatic 1/10 dilution of the sample. Bypassing manual sample dilutions reduces dilution errors and cross-contamination, making DiluCell™ ideal for GLP. Combined with small sample volume requirements and bubble-free filling, the DiluCell™ allows for convenient spectrophotometric analysis from 230 - 950 nm (UV-compatible). This feature enables the user to incorporate the DiluCell™ cuvettes for all standard applications in addition to their UV-Vis-based measurements.

Operating instructions can be found in the manual under “*DiluCell™ technology*”.

The DiluCell™ 10 can be ordered at Implen under the article number: DC 10 UV-Vis.

Connectivity

▪ Wi-Fi® (Wifi)

The Implen OD600® is equipped with Wi-Fi®, which can be used as a Wi-Fi® network or as a Wi-Fi® Hotspot. The Wi-Fi® can be used to integrate into a LIMS system via a REST API.

To connect to the LIMS, please contact Implen Support.

To connect to the device's hotspot or to connect your device to Wi-Fi®, please refer to the "*Network*" section in the manual.

▪ USB-A

There is a USB-A port on the front and rear panel of the Implen OD600® which is compatible with standard portable USB 2.0 storage devices for direct data transfer in Excel (.xls) and .csv formats. Encrypted USB flash drives are not compatible.

Note: We recommend using FAT/FAT32 formatted 2.0 USB flash drives.

Note: If the USB stick is not working, please re-format it to FAT32 using a Windows PC.

▪ LAN/Ethernet

There is an Ethernet (LAN) connection port on the rear panel of the instrument which enables the Implen OD600® to integrate with a LIMS system via a REST API.

To connect to the LIMS, please contact Implen Support.

LAN functionality is activated automatically when a LAN cable is connected.

Note: If the LAN-cable was not detected by the device restart the device and plug the calbe in again.

Note: The maximum LAN cable length is 10 meters. Bit rate is 100 Mbit/s.

Implen OD600® Specifications

Specification	Details	
Wavelength	600 nm	
Bandwidth	18 nm	
Photometric Range	OD600: 0 - 4 A	McFarland: 0 - 16 MFU
Accuracy	OD600: @ 1A: < ±0.01 A	McFarland: @ 0 - 8 MFU: ≤ ±0.1 MFU
Reproducibility	OD600: @ 1A: < ±0.002 A	McFarland: @ 0 - 8 MFU: ≤ ±0.05 MFU
Cuvettes/Tubes	DiluCell™, 10mm macro and semi-micro cuvettes; 10, 12, 16, 18 mm glass/falcon tubes - without the need for adapters. Wipeable sample compartment	
Optical height	8.5 mm	
Displayed information	Sample list function: View sample data onboard incl. Abs value, MFU value, sample name, time and date of measurement	
Display size	3,5" glove compatible touchscreen, 480 x 320 pixels	
Memory	4 GB internal memory, access to all previously measured sample data in sample list feature	
Time and date data	Manual or automatic time setting	
Connection capability/output	USB-A port for USB stick data export Wi-Fi® and Ethernet/LAN for LIMS integration; USB-C for Charging	
Battery life	Rechargeable, continuous use for ~8 hours	
Power requirements	External power adaptor and internal rechargeable battery	
Data format	.xls and .csv format exportable; For LIMS: JSON	

General Specifications

Light Source	LED
Detector	Photodiode
Size	140 mm x 110 mm x 70 mm
Weight	0.6 kg
Display	480 x 320 pixels; touchscreen glove compatible
Operating Voltage	100-240 V, 50/60 Hz, 20 W, 0,7 A
Certifications	CE, IEC/EN 61010-1:2010 and EN 61326-1:2021
Battery Pack	Rechargeable lithium ion battery; 3.7V, 4400mAh Operation time: up to 8 h; min. charging cycles: 500
Battery Certification	UN38.3 transport test

Features and specifications are subject to change without notice.

2. Getting Started

Implen OD600® Installation

Note: Translations of the Safety Information and Unpacking and Positioning are available in Appendix 1 (Czech), Appendix 2 (Dutch), Appendix 3 (French), Appendix 4 (German), Appendix 5 (Italian), Appendix 6 (Polish), Appendix 7 (Spanish), Appendix 8 (Swedish), Appendix 9 (Portuguese), Appendix 10 (Bulgarian), Appendix 11 (Croatian), Appendix 12 (Danish), Appendix 13 (Estonian), Appendix 14 (Finnish), Appendix 15 (Greek), Appendix 16 (Hungarian), Appendix 17 (Irish), Appendix 18 (Latvian), Appendix 19 (Lithuanian), Appendix 20 (Maltese), Appendix 21 (Romanian), Appendix 22 (Slovak), Appendix 23 (Slovenian), Appendix 24 (Chinese), and Appendix 25 (Japanese).

▪ Safety Information

Before commencing installation, please take time to familiarize yourself with warning labels and symbols on your instrument and their meaning. These are to inform you where potential danger exists or particular caution is required.

 Direct current

Do not open the instrument as this can expose the operator to electrical power and delicate optical setup or damage the instrument.

Do not use damaged power cords, accessories, and other peripherals with your Implen OD600®. Use only the delivered and specified power supply/charger.

Do not expose the Implen OD600® to strong magnetic, electrical fields, water or any type of liquid as heavy rain or moisture.

Do not put the instrument into fire, as it may swell or explode (battery). Do not store at or use near any type of heat source, especially temperatures above 50°C.

Do not leave your Implen OD600® on your lap or near any part of your body to prevent discomfort or injury from heat exposure.

Do not place objects on top of the Implen OD600®.

The Implen OD600® has to be switched off during transportation.

The symbol  on the product, or on the documents accompanying the product, indicates that this appliance may not be treated as household waste. Instead it shall be handed over to the applicable collection point for the recycling of electrical and electronic equipment. Disposal must be carried out in accordance with local environmental regulations for waste disposal.

▪ Unpacking and Positioning

Check the contents of the package against the delivery note. If any shortages are discovered, inform your supplier immediately.

Inspect the instrument for any signs of damage caused in transit. If any damage is discovered, inform your supplier immediately.

Ensure your proposed installation site conforms to the environmental conditions for safe operation: indoor use or dry environment.

The Implen OD600® is designed and developed for use in a controlled, secure environment. This is described as a laboratory environment with restricted physical access within a company. The Implen OD600® may only be accessed by qualified users and may only be operated within a local company network. If a secure environment is not ensured, the Implen OD600® does not meet the following requirements of the Radio Equipment Directive (RED DA): ACM, AUM, SUM, SSM, SCM, CCK, GEC.

Note: Do not expose your Implen OD600® to liquids, rain, moisture or dusty environments.

Operational temperature range: 15 - 40°C

Storage temperature range 5-50°C

Transport temperature range: -20°C to 50°C

If the instrument is subjected to extreme temperature changes, it may be necessary to allow the instrument to equilibrate and to turn the instrument off and then on again once thermal equilibrium has been established (~2 - 3 hours).

Maximum relative humidity (non-condensing) of 80% and up to 31°C decreasing linearly to 50% at 40°C.

The instrument must be placed on a stable, level surface that can support 0.6 kg and so that air can circulate freely around the instrument. Ensure while powered on that no materials reduce air circulation.

The equipment should be positioned such that in the event of an emergency the main plug can be easily located and removed.

Carry the instrument always by holding the main corpus of the instrument.

The equipment must be connected to power with the USB-C charger – output: 20W max., input: AC 100-240V, 50/60Hz, 0.7A. The power supply/cord supplied by Implen is a USB-C adapter with 5V. Please charge the battery for at least for 3 hours prior to the first use.

Turn the instrument on using the power button on the touch-panel after it has been plugged in. The instrument will perform a series of self-diagnostic checks.

Please read the complete user manual prior to use.

Please contact original supplier immediately if technical or sample handling difficulties are experienced.


Note: If this equipment is used in a manner not specified or in environmental condition not suitable for safe operation, the protection provided by the equipment may be impaired and the instrument warranty voided.

3. Implen OD600® Basics

The Implen OD600® product offers a complete solution for standard cuvette/test tube applications. Glass or plastic test tubes/falcon tubes (10, 12, 16, 18 mm diameter - no adapters are required) as well as all standard macro or semi-micro cuvettes can be used (quartz, glass, or plastic cuvettes with a center height of 8.5 mm).

Applications Overview

The Implen OD600® is equipped with pre-programmed applications, the OD600 and McFarland methods. To select an application, tap on the corresponding icon and the method will open immediately.

Method Icons

Description

OD600


Measures cell density at ~600 nm

McFarland

Estimate cell concentration by comparing turbidity to standard solutions

Icons

Icon	Name	Action
	Battery Charging	Battery power plugged and full
	Battery Full	Battery charged 100 %
	Battery Medium Full	Shows the actual battery status
	Battery low	Battery critically low, please re-charge immediately
	Wi-Fi® Network	Wi-Fi® network active
	Wi-Fi® Hotspot	Wi-Fi® Hotspot active
	Wi-Fi® Network & Wi-Fi® Hotspot	Wi-Fi® Network & Wi-Fi® Hotspot active
	LAN	LAN connection active
	Menu	Opens the Menu page
	Settings	Opens the Settings page
	Home Screen	Returns to the home screen with application icons for method selection.
	Back	Returns to the previous page
	Delete Data	Opens a delete dialog pop up
	Tube Selection	Opens a menu of container types from which you can choose your desired container type
	Sample List	Opens a sample list of sets containing all measurement data with further action options including delete and export of data files.

Buttons

BLANK

Prior to any sample measurement at the start of a new method a blank measurement of either water or the buffer of the samples is required to give the Implen OD600® a reference of what zero should be. It is recommended to re-apply the blank and measure it as a sample to ensure the graph of the blank spectra is a flat line. To ensure accurate measurements, we recommend blanking the device before starting a measurement if 30 minutes or more have passed since the last blanking.

SAMPLE

To initiate the measurement of the sample press the sample button. The data will be stored automatically in the sample list function.

General Features

Sample List

This section will guide you through using the Sample List to view, manage, and export your sample measurement data. All data will be auto-saved to the Sample list.

1. Accessing the Sample List:

From the main OD600 or McFarland Application pages, tap on the Sample List icon.

A new page will open, displaying a list of your previously measured samples.

The Sample List is scrollable for easy navigation.


Use arrows next to the table to move through the list.

2. Sample Measurement Sets:


Sets are listed in reverse chronological order, labeled with a Set number and creation date/time.

Measurements are grouped into sets of experiments.

Each set represents a session of using the app, from opening to closing of the application.

To select a sample or set, tap on the corresponding row. Selected items will be highlighted in blue.

When selecting a set, all samples within that set are automatically included.

Four samples/sets can be selected at once. Selected samples can be deleted or exported for further analysis. Use the ALL button to quickly manage selections.

If your device enters deep sleep mode, a new set will be automatically started.

Note: If the data is selected manually, be aware that only the four samples displayed on the screen will be selected. The selection will be lost after scrolling down to the next samples.

3. Exporting Results:

Please insert a USB storage device to proceed with the export.

Select sets/samples to export, then press the EXPORT button.

You can export your data in the following formats: .csv or .xls (default). Once a format is selected, tap the EXPORT button to initiate the export process.

The data will be automatically saved in a new folder created on the USB stick.

Note: When exporting data from the sample list, the Excel format allows for exporting fewer than 200 samples at a time. For larger datasets, the CSV format should be used, as it does not impose any limit on the number of samples that can be exported in a single operation.

Note: If the USB stick is not working, please re-format it to FAT32 using a Windows PC.

Sample Handling Tips

The cuvette holder is compatible with glass or plastic test tubes/falcon tubes (10, 12, 16, 18 mm diameter - no adapters are required) as well as all standard macro or semi-micro cuvettes (quartz, glass and plastic cuvettes with an optical height of 8.5 mm).

It is also possible to use cuvettes with 5 mm, 2 mm or 1 mm path lengths, but there may be an adapter necessary. Please ask your cuvette supplier for suitable adapter.

The minimum volume for accurate measurements depends on the cuvette/tube type used; it is necessary for the light path to pass through the sample for accurate measurements.

Note: Do not overfill the cuvette for risk of spillage.

Note: The cuvette holder is not removable.

Note: Do not use aggressive solvents such as strong acids, bases, or organic solvents. Highly concentrated acids and bases are also not recommended. After each measurement, wipe the sample surface with a lint-free laboratory wipe. For detailed information on compatibility of specific reagents or solvents, please contact Implen support at support@implen.de.

Operating Conditions

When operated under anaerobic (oxygen-free) conditions, the presence of carbon dioxide (CO₂) can cause long-term corrosion of metallic components, particularly carbon or low-alloy steels.

To minimize corrosion risk:

- Use nitrogen (N₂) as the inert cover gas wherever possible.
- Ensure the gas is dry (low humidity) and the system is free of condensate.
- Avoid prolonged exposure to CO₂-containing or acid-forming gases.
- Maintain neutral to mildly alkaline conditions if a liquid phase is present

Power Management Features

Battery Operation

A full battery pack can be operated for approximately 8 hours within normal usage. When the battery is close to empty, a single warning popup will appear 1 minute before shutdown. At this point, it is essential to connect the device to a power source immediately.

The battery pack is designed for at least 500 full charging cycles. Afterwards the capacity and operation time of the battery may vary. A fully charged battery pack will self-discharge (switched off device) over a time of 10 – 14 days. After this time recharging of the battery pack is necessary.

Note: If the battery status has not reached 100% after the maximum charging time, keep the device turned on and connect it to the main power. The battery status should then update.

Auto-Shutdown Features

If the device is not charged and remains idle on a desk without any operation, it will follow a specific power-saving sequence:

- After 10 minutes: The screen will automatically dim to conserve energy.
- After 29.5 minutes: A popup notification will appear, stating that the device is about to shut down.
- After 30 minutes: The device will shut down automatically if no operation has occurred.

However, if the device is connected to a charger while not in use, it will remain powered on indefinitely, without initiating the shutdown sequence.

Note: If the Implen OD600® shuts down during a measurement; the data of this measurement are lost.

4. Implen OD600® Applications

The Implen OD600® comes with two pre-programmed applications, the OD600 application and the McFarland application. Any application method can be selected by tapping the icon once.

OD600 Application

Method Overview

The growth of bacteria in liquid culture media is commonly monitored by measuring the optical density at 600 nm (OD600) in small samples taken from the cultures. OD600 measurements are typically used to determine the stage of growth of the bacterial culture, thereby ensuring that cells are harvested at an optimum point that corresponds to an appropriate density of live cells. Growth of bacterial cells typically progresses through a series of consecutive phases including: lag, log, stationary and decline (see figure 1). In general, cells should be harvested towards the end of the log phase using the optical density of the samples to determine when this point has been reached. Since optical density in the case of OD600 measurements results from light scattering rather than light absorption, this value varies depending on the type of bacterial cells in the culture in terms of size and shape. Cells are routinely grown until the absorbance at 600 nm (known as OD600) reaches approximately 0.4 prior to induction or harvesting. A linear relationship exists between cell number (density) and OD600 up to an absorbance value of 0.6, approximately.

As mentioned above, for turbid samples such as cell cultures, the absorbance measured is due to light scattering, and not the result of molecular absorption. Since the extent of scattering is affected by the optics of the system (distance between the cell holder and instrument exit slit, monochromator optics, slit geometry, etc.), different spectrophotometer types will tend to give different OD600 readings for the same turbid sample. Therefore, if results from different spectrophotometers are to be compared, they must be normalized first using appropriate calibration curves.

A calibration curve can be constructed by comparing measured OD600 to expected OD600. Expected OD600 is determined by counting cell number using an alternative technique (for example microscope slide method) and converting to OD600 using the rule of thumb that $1 \text{ OD600} = 5 \times 10^8 \text{ cells/ml}$ for *E. coli*.

Note: The use of 10 mm path length disposable cells is recommended for optical density measurements of cell culture solutions. To prevent the suspension settling too quickly and giving an OD reading that changes with time, glycerol should be added to the sample.


Figure 1: Bacterial growth curve

Measurement Protocol

1. Getting started

To begin using the OD600 application, simply tap on the icon displayed on your device's Home screen.

2. Perform Blank Measurement


BLANK

Add a blank solution to a cuvette and ensure that the filling volume (8.5 mm optical height) is sufficient to allow light path to pass through the solution.


Starting a new measurement will create a new set of data that is automatically saved to the sample list.

Note: The sample button is inactivated until a blank has been measured.

Note: The light path is indicated with two white arrows on the device. Make sure that the cuvette is inserted into the device in the right direction. The arrows on the cuvette should match the arrows on the device.

Initiate a blank measurement with the blank button.

Upon completion of the measurement, remove the cuvette.

Note: To ensure accurate measurements, we recommend blanking the device before starting a measurement if 30 minutes or more have passed since the last blank.


3. Perform a Sample Measurement

Add your sample solution to a new cuvette and ensure that the filling volume (8.5 mm optical height) is sufficient to allow light path to pass through the solution.

Insert cuvette/tube into the cell holder.

Initiate a sample measurement with the sample button. Upon completion of measurement remove the cuvette.

Apply further samples.

Note: Handling Out-of-Range Measurements - If a sample measurement is out of range (>4A), a warning popup will appear, instructing you to repeat the measurement and dilute the sample accordingly.

Note: If the light path is blocked by an obstacle, the error code Err4 will appear - Please remove the obstacle from the light path to proceed with the measurement.

4. Saving and Accessing Measurements

All sample measurements and sets are automatically saved and can be accessed through the Sample List feature in the respective app. Data can be exported through the sample list function using a USB-Stick.

Calculations

$$OD600 = A_{600} * D$$

OD600	Optical density at 600 nm
A ₆₀₀	Absorbance at 600 nm (10 mm path)
D	Dilution factor

$$Cells/ml = A_{600} * D * cf * multiplier$$

A ₆₀₀	Absorbance at 600 nm (10 mm path)
D	Dilution factor
cf	Correction factor for spectrophotometer
multiplier	Multiplier of sample

McFarland


Method Overview

The McFarland standard is a technique used to estimate the density of bacterial suspensions by visual comparison with a series of turbidity standards. This method is particularly useful when direct cell counting is not feasible or practical. The McFarland standard is commonly employed in microbiology laboratories for various applications such as antimicrobial susceptibility testing and bacterial inoculum preparation.

The procedure involves preparing a series of standard turbidity solutions known as McFarland standards. These standards have specific turbidity levels corresponding to different optical densities. The most common McFarland standard is 0.5, which is roughly equivalent to 1.5×10^8 colony-forming units (CFU) per milliliter for most bacterial species.

A sample suspension is compared to the standard curve established using McFarland standard solutions with known turbidity values. This provides an estimate of the cell density in the suspension.

A calibration curve for the 16mm test tube is saved on the device as a standard and can be used as the default setting.

Note: The default 16mm tube calibration curve was generated using the McFarland standards set from Company ProLab, Serial Number 395123, traceable to 260-140 SD2350.


Note: Please be aware that standards can vary between batches and supplier. For optimal accuracy, it is recommended to create your own calibration curve using your own standards, as these can differ and may expire over time.

Container selection & Calibration

This guide provides step-by-step instructions for selecting tubes and calibrating your instrument for accurate measurement of turbidity using McFarland units (MFU). Proper tube selection and calibration are essential for obtaining reliable results.


1. Container Selection

Before calibrating your instrument, ensure the correct tube/cuvette type is selected. Press the button indicating the tube that is currently in use.

Note: The tube type selection in the McFarland app is intended solely for labeling purposes and does not dictate the specific tube size required for calibration with the standard curve. It is not mandatory to use the indicated tube type; however, for accurate and reliable results, it is essential to use the same tube type consistently across calibration, blank, and sample measurements.

A menu table will appear with options for tube types, including 10 mm cuvette, 10 mm test tube, 12 mm test tube, 16 mm test tube (default), and 18 mm test tube.

The default tube type for McFarland measurements is the 16 mm test tube. Implen has factory-set a 16mm test tube calibration curve, which can be overwritten, which will be reset after a factory calibration.

The blue-highlighted tube indicates the one currently selected calibration. The tick-box indicates the selected calibration, which will be memorized by the device.

Tap on the corresponding tube type to select a new tube-setting.

If no calibration has been done for the selected tube type, a popup will state "No calibration data available. Start new calibration?"

You can then choose to exit or proceed with calibration.

Note: The default 16mm tube calibration curve was generated using the McFarland standard set from Company ProLab, Serial Number 395123, traceable to 260-140 SD2350.

Note: Please be aware that standards can vary between batches and supplier. For optimal accuracy, it is recommended to create your own calibration curve using your own standards, as these can differ and may expire over time.

2. Calibration Process

After selecting the tube type, a new page will open displaying the date of the previous calibration (if any).

Note: If a tube type with an existing calibration is selected, be aware that a new calibration will override the existing calibration!

Two actions are available: SELECT and CALIBRATE.

Note: If there is an existing calibration available, the options will be: CALIBRATE and EXIT.

SELECT leads you to the measurement window.


If CALIBRATE is selected, a calibration page is opened. Follow the instructions on the calibration page to perform blank and standard measurements:

Measure the blank solution.

Ideal Case: The 0.00 McFarland standards should be included in the calibration set. It represents the baseline or zero turbidity and helps ensure that the instrument is properly calibrated to detect small differences in turbidity.

Alternative for 0.00 Standard: If a 0.00 McFarland standard is not available, prepare a water blank in a cuvette of the same diameter as those used for measurements. This water blank serves as a substitute for the 0.00 standard and helps in accurately setting the baseline for measurements.

Click on Cal Curve to start with the measurements to create a new standard curve.


Measure the standard solutions one after the other in ascending order of turbidity i.e. 0.5, 1, 2, 3, and 4.

If you do not have a full set of 5 standards, you can skip some. A minimum of 2 Standards is required.

Note: A minimum of two standards is necessary to create a new calibration curve.

After measuring all standards, the calibration curve will appear. After all measurements are completed, the calibration data will be automatically saved. The R2 value of the curve will be displayed. Click "RUN SAMPLE" to proceed with your measurements.

Note: Ensure all cuvettes and test tubes are clean and free from any residues or scratches that might affect the measurements.

3. Troubleshooting

When troubleshooting calibration issues, please remember that McFarland standards can vary between batches and suppliers. This variability can affect the accuracy of your measurements. For best results, it may be necessary to gently invert the standards before use to ensure they are well-mixed. Additionally, consider creating a new calibration curve using your own standards, as they may differ from the default curve and could expire over time.

If you encounter any other issues during tube selection or calibration, refer to the troubleshooting section of the instrument's manual or contact technical support for assistance.

Measurement Protocol

1. Getting started

To begin using the McFarland application, simply tap on the icon displayed on your device's Home screen.

2. Select a tube type

The default calibration curve provided is for the 16 mm test tube.

If you wish to use a different calibration curve or test tube, click on the button displaying the tube type and select a different one (indicated by yellow arrow).

If no custom calibration curve has been saved previously, you will need to establish a new calibration curve. For instructions on this process, please refer to the manual section *Container Selection & Calibration* above.

Note: The default 16mm tube calibration curve was generated using the McFarland standard set from Company ProLab, Serial Number 395123, traceable to 260-140 SD2350.

Note: Please be aware that standards can vary between batches and supplier. For optimal accuracy, it is recommended to create your own calibration curve using your own standards, as these can differ and may expire over time.


3. Perform Blank Measurement


Add a blank solution to a cuvette and ensure that the filling volume (8.5 mm optical height) is sufficient to allow light path to pass through the solution

Note: The light path is indicated with two white arrows on the device. Make sure that the cuvette is inserted into the device in the right direction. The arrows on the cuvette should match the arrows on the device.

Note: The sample button is inactivated until a blank has been measured.

Initiate a blank measurement with the blank button.


Upon completion of the measurement, remove the cuvette.

4. Perform a Sample Measurement

Add your sample solution to a new cuvette and ensure that the filling volume (8.5 mm optical height) is sufficient to allow light path to pass through the solution.

Insert cuvette/tube into the cell holder.

Initiate a sample measurement with the sample button. Upon completion of measurement remove the cuvette.

Note: Handling Out-of-Range Measurements - If a sample measurement is out of range (too high absorbance or lightpath blocked), a warning popup will appear, instructing you to repeat the measurement and dilute the sample accordingly.

Note: If the light path is blocked by an obstacle, the error code Err4 will appear - Please remove the obstacle from the light path to proceed with the measurement.

5. Saving and Accessing Measurements

All sample measurements and sets are automatically saved and can be accessed through the Sample List feature in the main menu. Data can be exported through the sample list function using a USB-Stick.

Calculations

The McFarland concentration is calculated using a basic linear regression standard curve.

This plot is described by the equation of a line: $y=mx+b$ where:

- y is the response variable (measured response),
- x is the explanatory variable (concentration),
- m is the slope of the line,
- b is the y-intercept.

A line of best-fit is then used for concentration determination. Once the equation/curve is programmed, the unknown concentration/MFU value is calculated based on these values.

DiluCell™ Technology

Method Overview


The DiluCell™ is a disposable cuvette designed for virtual dilution of cuvette-based samples, utilizing the principle of the Lambert-Beer law, which empirically relates the absorption of light to the properties of the sample.

Lambert-Beer Law:

$$Abs = \log \frac{I_0}{I} = \epsilon \times c \times d$$

Where:

I_0	=	incoming light
I	=	outgoing light
ϵ	=	extinction coefficient
c	=	concentration
d	=	pathlength (1mm @ 8.5mm center height)


Built-in Virtual Dilution

Traditional cuvettes require manual sample dilution and volumes of at least 1 ml.

The innovative design of the DiluCell™ reduces the light path from the standard 10mm to just 1mm, thereby enabling a **built-in 1:10 virtual dilution** in accordance with the Lambert-Beer Law.

This allows you to work with smaller sample volumes (only 300 µl), while eliminating manual dilution steps, saving time, and minimizing the risk of dilution/pipetting errors and cross-contamination.

Operating Instructions

1. Correct Insertion

Place the DiluCell™ in the cuvette holder with the wide window facing the light beam. The Implen OD600® features a directional arrow to assist with correct alignment.

2. Sample Loading

The geometry of DiluCell™ supports bubble-free filling. After loading, check visually for air bubbles. If bubbles are present, tap gently to eliminate them.

3. Volume Requirement

Use 300 µl of sample for optimal performance.

4. Data Handling

Simply multiply your absorbance values by 10 post-measurement or set your instrument to apply a 10× dilution factor.

▪ Important Usage Notes

Note: DiluCell™ will not work in any other orientation; do not mix other cuvette types within a DiluCell™ measurement series (for blank, etc.).

Note: Single Use Only. It is not recommended to reuse DiluCell™ due to: high probability of cross contamination and high probability of air bubbles when refilling the cell.

Note: Do not autoclave the DiluCell™.


Note: For solvent compatibility please contact your local Implen partner or the Implen Team at www.implen.de.

Specifications

Parameter	DiluCell™ 10
Wavelength Range	340 nm – 950 nm
Optical Pathlength	1 mm (at 8.5 mm center height)
Volume Requirement	300 µl
Tolerance	±10%
Dilution Factor	1:10
Product Code	DC 10
Packaging	96 disposable cuvettes per pack

5. Preferences


System preferences can be set by selecting preferences (⚙️) on the home screen. The Main Menu includes: Settings, Storage, About and Help.

Main Menu Options

Settings

Upon selecting the "Settings" option from the preferences menu, a submenu will appear with several configuration options. These options include Date and Time, Network, Display, Language, and System setup.


▪ **Date and Time**

Within the Date and Time settings, you can choose between manual and automatic time setting. Manual time setting allows you to adjust the date and time manually using buttons on the screen. Automatic time setting synchronizes the device's time with an external time server when connected to the internet.

A timezone can be selected through the Implen OD600® webpage (see instructions under “Network > Implen OD600®” webpage below)

Note: If your device is not connected to a network, the automatic time zone cannot be activated. To activate the automatic time zone, please connect the device to Wi-Fi® or LAN and change timesetup through your internet browser on the setup webpage.


Note: If you want to change the time zone through the webbrowser using the device-hotspot, ensure Wi-Fi® and Wi-Fi® Hotspot are activated before pressing save, otherwise the device will disconnect from the computer before changing the time zone.

Note: If the time zone or date/time is adjusted through the web browser via the Implen OD600 webpage, the changes may take several seconds to take effect, and the updated date/time on the device will be reflected within the next minute.

▪ **Network**

The Network submenu provides options for configuring network settings such as Wi-Fi®, Wi-Fi® Hotspot, and LAN.


Wi-Fi® Hotspot

Follow these steps to connect your device to the Wi-Fi® Hotspot and configure the Wi-Fi® settings using the Implen OD600® webpage.

▪ **Step 1: Activate Wi-Fi® Hotspot on Implen OD600®**

1. Navigate to the Network Menu on the Implen OD600®:

Go to **Main Menu > Settings > Network > Wi-Fi® Hotspot.**

1. Activate the Wi-Fi® Hotspot.

The SSID and IPv4 Address will be displayed in the Wi-Fi® Hotspot menu.

▪ **Step 2: Connect Your Device to the Wi-Fi® Hotspot**

For Smartphones:

1. Scan the QR-code located on the bottom of the Implen OD600® to connect your smartphone to the device's hotspot.
2. Open a web browser on your smartphone and enter the IPv4 Address: 192.168.4.1.
3. This will open the device's network settings interface.

• **For Computers:**

1. Connect your computer to the Implen OD600® Wi-Fi® Hotspot.
2. Use the password **“Implenuser”** to connect to the hotspot.
3. Open a web browser on your computer and enter the IPv4 Address: e.g. 192.168.4.1.
4. This will open the Implen OD600® webpage to change the device's network settings interface

▪ **Step 3: Configure Network Settings through the Implen OD600® webpage**

Configuration of network settings via the Implen OD600® webpage is explained in the section “Implen OD600® webpage” below.

Wi-Fi® (Wifi)

▪ **Connecting Implen OD600® to Wi-Fi® for the First Time**

1. **Activate the Wi-Fi® Hotspot on the Implen OD600®:**

- Navigate to **Main Menu > Settings > Network > Wi-Fi® Hotspot** and activate the hotspot. Refer to **Step 1 in the Wi-Fi® Hotspot Setup** section.
- The device will display the SSID and IPv4 Address in the Wi-Fi® Hotspot menu.

2. **Connect Your Device to the Wi-Fi® Hotspot:**


- **For Smartphones:** Scan the QR-code on the bottom of the Implen OD600® device to connect, then open a web browser and enter the IPv4 Address

- **For Computers:** Connect to the hotspot using the password “**Implenuser**”, then enter the IPv4 Address in a web browser on your computer. This will open the Implen OD600® webpage.

3. **Configure Wi-Fi® Network Settings via the Implen OD600® Webpage:**

- Once connected to the hotspot, the Implen OD600® webpage will open.
- Enter the desired Wi-Fi® network’s SSID (Wi-Fi® Name) and password.
- Click **Save** to apply the settings. Refer to the section “Implen OD600® webpage” below for more information.
- Check the configuration on your OD600 device

▪ **Activating Wi-Fi® on Implen OD600® After Initial Setup**

Once the Implen OD600® has been connected to a Wi-Fi® network, you can activate the Wi-Fi® function directly without reconnecting through the hotspot:

1. **Navigate to the Wi-Fi® Menu:**

- Go to **Main Menu > Settings > Network > Wi-Fi®**.

2. **Activate Wi-Fi®:**

- Toggle the Wi-Fi® function to "On."
- The device will automatically connect to the previously configured Wi-Fi® network.

Note: Ensure that only the Wi-Fi® or LAN function is active at any given time to avoid network conflicts. Refer to the Network Overview section for additional guidance.


LAN/Ethernet

LAN functionality is activated automatically when a LAN cable is connected.

The Implen OD600® Webpage can be accessed through entering the IPv4 Address.


Note: If the LAN-cable was not detected by the device, restart the device and plug the cable in again.

Note: Be cautious not to activate multiple network functions simultaneously to avoid errors.

Implen OD600® Webpage

The Implen OD600® webpage allows you to configure device settings, including date & time, Wi-Fi®, and Wi-Fi® hotspot configurations. Additionally, you can view the End User License Agreement (EULA) and licenses.


This webpage can be accessed via Wi-Fi® or LAN connections. Follow the steps below to connect and use the Implen OD600® webpage.

Accessing the Implen OD600® Webpage

▪ **Option 1: Connecting via Wi-Fi® Hotspot**

1. **Activate Wi-Fi® Hotspot on Implen OD600®:**

- Navigate to the device's Network Menu:

Go to **Main Menu > Settings > Network > Wi-Fi® Hotspot.**

- Ensure the Wi-Fi® Hotspot is activated.
- The SSID and IPv4 Address will be displayed in the Wi-Fi® Hotspot menu.

2. **Connect a Smartphone:**

- Scan the QR-code located on the bottom of the Implen OD600® to connect your smartphone to the device's hotspot.
- Open a web browser on your smartphone and enter the IPv4 Address: 192.168.4.1.
- This will open the device's network settings interface.


3. **Connect a Computer:**

- Use your computer to connect to the Implen OD600® Wi-Fi® Hotspot.
- Use the password **“Implenuser”** to connect to the hotspot.
- Open a web browser on your computer and enter the IPv4 Address: 192.168.4.1.
- This will open the device's network settings interface.

4. Configure Network Settings:

- In the network settings interface:
 - Enter the Network SSID (Wi-Fi® Name) and password of the desired Wi-Fi® network to connect the Implen OD600®.
 - The Wifi and Wifi Hotspot functions can be activated in this menu.
 - Click **Save** to apply changes.

5. Adjust Date&Time Settings:

- Use the interface to adjust the device's date & time settings. Chose the desired timezone and click **Save**.
- **Option 2: Connecting via LAN**
 1. **Connect via LAN:**
 - Connect the Implen OD600® to LAN.
 - On your device, navigate to the LAN menu and obtain the IPv4 Address.
 - Open a web browser on your computer and enter the IPv4 Address provided in the LAN menu.
 - This will open the device's network settings interface, allowing you to configure settings as needed.

▪ **Display**


In the Display settings, you can adjust the brightness levels of the device's screen using a slider bar.

The slider bar offers smooth transitions with 10-20 brightness increments.

The default brightness is 50%.

Any changes to the brightness will be saved and set as the new default when the device is turned on and off.


After performing a factory reset, the brightness will revert to the original default of 50%.

▪ **Language**


Language settings allow you to select the preferred language from the available options.

After selecting a language, a popup will prompt you to restart the device for the changes to take effect.

▪ System Setup

The System Setup menu includes options for software updates, creating log files, factory reset, and data backup.

Software Update:

Please insert a USB-Stick into the device, containing a software update from the Implen homepage. Software updates can be found here: <https://www.implen.de/downloads/>

Create a log file:

Please insert a USB-Stick into the device. The Log file will be exported to the USB-Stick after confirmation.

Factory Reset:

After the factory reset, all self-acquired data will be lost.

The device will restart and show the license agreement page again.

The Wi-Fi® function of the device will be turned off. To turn it back on, use the Wi-Fi® Hotspot (see “Network”).

If the device is not connected to the Wi-Fi® anymore, time and date need to be adjusted manually.

Backup all data:


A backup of all data can be done via USB export.

Storage

Access the Storage option from the main menu to view available storage space and usage.


You can clear all measurement data from the device's memory using the Clear All Measurement Data button.

Confirmation popups will appear before clearing data, and warnings will be displayed for low or full storage.

About


The About section provides device information such as the device name, serial number, and firmware version.

Help

Access helpful resources such as the user manual and contact details for support.

QR codes are provided for accessing the user manual and contacting support.

6. Troubleshooting

Error Indication	Solution
A negative reading is obtained.	In normal measurements the test sample has a positive Absorbance compared to that of the Reference. Negative readings will be obtained if the Reference and Test cuvettes are mixed up.
The OD600 value is different to that obtained on another instrument in the lab.	When you measure turbid solutions you do not measure the absorbance/transmittance of light at the detector, you measure the amount of scattered light that reaches the detector. Thus optical geometry is very important - the further the distance from the sample to the detector, the greater the effect of the scattered light. Thus instead of harvesting at 0.4 OD, for example, you have to do it at 0.8 OD. A simple conversion factor can be calculated from the OD 600 of your existing instrument compared to that of the Implen OD600 instrument. For further information see also Technical Note #8 OD600. Note: the technical notes can be found on implen.de
Device shuts down and restarts after ESD-impulse	If the device shuts down after an ESD-impulse, it will restart and move back into the home screen.
Device shuts down and restarts: Light Path Obstruction Error	If the device crashes and a sample is still inside, the light path for the self-test is blocked. Solution: Remove the sample to clear the light path, then blank again.
The USB Stick is not detected by the device	If the USB stick is not working, please re-format it to FAT32 using a Windows PC.
The selected data was not exported	In the sample list screen, only the four samples in the visible screen can be selected at a time. If you scroll to the next page, the samples from the previous page will be unselected. If you want to export all data, click on <i>ALL</i> .

Note: the technical notes can be found on implen.de

Self-test

The Implen OD600® Self-test is performed automatically every time the instrument is powered on. If the instrument passes the self-test, the home screen will appear after a few seconds.

If the instrument does not pass the test, a window with the message explaining the reason for the failed test along with the recommended solution will be shown. If the OK button is selected, the window is closed and the home screen will be shown with an error message at the top. The error message is also shown on top of all other screen and all performed measurements including the saved files of those measurements. If the self-calibration test fails please contact the Implen Support Team.

7. Important Error Messages

Error	Information	Action
Err1	Blank absorbance to low	Blank absorbance to low. Please contact Implen Support. Err1
Err2	Blank absorbance to high	Blank absorbance to high. Remove the obstacle in the light path! Otherwise contact Implen Support. - Err2
Err3	Sample absorbance to low	Sample absorbance to low. Please contact Implen Support. - Err3
Err4	Sample absorbance to high	Sample absorbance to high. Remove the obstacle in the light path! Otherwise contact Implen Support. - err4
Err5	Hardware Fail	Hardware failed. Contact Implen Support. - Err5
Err6	General Measurement Error	Generic error! Please contact Implen Support. - Err6

8. Important Warning messages

Warning Message	Solutions
No USB device detected. Please connect a USB flash drive to the USB port.	Disconnect and reconnect the USB device. Use a FAT/FAT32 formatted USB 2.0 flash drive. If the USB stick still doesn't work, reformat it to FAT32 using a Windows PC.
Battery cannot be charged. - [number]	Ensure the device is within the required operating temperature. If the issue persists, contact Implen support.
Battery critically low. Please connect the charger.	Connect the device to the charger immediately.
Calibration failed. A minimum of two standard points are required to create calibration curve.	Ensure that standards are significantly different from the blank absorbance value. At least two valid standards are needed to create a calibration curve.
No network connection. Automatic date and time synchronization not possible. Please connect to Wi-Fi® or LAN.	Connect the device to the internet via Wi-Fi® or LAN for automatic date and time synchronization. You can adjust the timezone through the Implen OD600 webpage under the "Network" section.
Generic error! Please contact Implen Support. - [number]	Contact Implen Support for assistance. Contact details can be found in the Main Menu under "Help."
Self-test failed, counts too high. Please contact Implen Support. - [number]	Restart the instrument and ensure the light path is clear. If the issue persists, contact Implen Support.
Self-test failed, counts too low. Remove the obstacle in the light path! Otherwise contact Implen Support. - [number]	Restart the instrument and ensure the light path is clear. If the issue persists, contact Implen Support.
Storage space full. No further measurements possible. Move data to an external storage and delete onboard results to free up space.	Export important data to a USB stick and delete existing data on the device. For full space recovery, perform a factory reset, but note that this will erase all existing data.
Low storage alert: Device performance may be affected. Move data to an external storage and delete onboard results to free up space.	Export important data to a USB stick and delete existing data on the device. If all data is secured, a factory reset will clear space but will delete all existing data.
Failed to write the data to the memory	Try writing the data again. If the issue persists, contact Implen Support.
Power touch timeout error. - [number]	Ensure the device is in the correct working environment. Restart the device. If the issue persists, contact Implen Support.

System Freeze

In case of significant Electrostatic Discharge (ESD) impulses it could happen that the device briefly shuts down and switches on again. The system will recover in the main home screen. However, if there is no reaction, please press the touch power key at the cover glass of the display and restart the system.

9. Maintenance

Maintenance-free technology

The Implen OD600® technology is maintenance-free. Regular maintenance and calibration is not necessary.

For facilities that are working according to national as well as international guidelines and standards including: Good Laboratory Practice (GLP), Good Manufacturing Practice (GMP), or ISO9000-9004; the proper performance of the device has to be tested and proven on a regular basis with individually set intervals. A certified Neutral Density filter can be used as an optional accessory. These filters are suitable for the control and documentation of the photometric accuracy of your system. Please contact your local Implen office or an authorized Implen partner for further information.

Support agreements that help to fulfill the demands of regulatory guidelines concerning GLP/GMP include: calibration certification using filters traceable to international standards (during production and quality control), certified engineers and calibrated test equipment, approved to ISO 9001 standard, automatic self-diagnostic calibration test during start of the Implen OD600®, result is documented in each data output file, and users are given the possibility to save a Implen Document Source (IDS) file (no data manipulation possible).

Software update

There is a button for software update. It is needed to plug in a USB flash drive with the software installation file to the USB port. The software update files are available on the Implen homepage: <https://www.implen.de/downloads/>.

Update procedure:

1. Copy the Software Update file to a USB flash drive
2. Plug in the USB flash drive to the USB port of the Implen OD600® device
3. Go to Settings > Setup > Software Update
4. Press the Update button. The update will start and once finished automatically reboot the device.

Note: Make sure that the OD600 is connected to power and that the power connection is not interrupted during the update.

Note: Make sure you unpack the ZIP file of the Software-download prior to updating your device.

Contact Implen

There is an option to contact Implen for instruments connected to the internet. For any support issues and questions please contact the Implen team directly:

Europe, Asia, South Pacific, Middle East, Africa

Phone: +49 89 72637180
Fax: +49 89 726371854
Email: support@implen.de
Website: www.implen.de

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Replacement Parts

▪ LED replacement

The Implen OD600® is equipped with a LED and should not need replacement until several years of use. In the unlikely event the lamp does need to be replaced, this should only be done by the manufacturer.

▪ Battery replacement

The battery can only be replaced by the manufacturer.

▪ Touchscreen/glass plate replacement

The touchscreen/glass plate can only be assembled or replaced by the manufacturer.

Cleaning and General Care

Switch off the Implen OD600® and disconnect the power cord prior to external cleaning. Use a soft wet cloth or dry microfiber cloth to clean all external surfaces. A mild liquid detergent may be used to remove stubborn marks.

The instrument requires little maintenance, but the following are considered good practice:

1. Keep the instrument clean and dry. Wipe off any spilt liquids immediately. Clean with a slightly damp cloth; a non-abrasive water-based soap or detergent may be used. The instrument may be wiped.
2. Remove the cuvettes from the instrument when not in use.
3. Store in a cool place away from corrosive chemicals or fumes.

Approved disinfectant solutions include: Apesin disinfection spray (Tana Chemi GmbH), Dr. Becher disinfection spray (Dr. Becher GmbH), and Lysoformin Spezial (Lysoform Dr. Hans Roseman GmbH).

Optional: The cuvette shaft can be cleaned using a lightly damp cloth or a cotton swab.

Note: Observe all necessary precautions if dealing with hazardous samples or solvents.

De-contamination Procedure

To decontaminate we recommend that the instrument is wiped with ethanol or other antibacterial detergent as required. The instrument can be sterilized using formaldehyde or ethylene oxide, but not with UV light (due to plastic degradation). The cuvette shaft can be cleaned using a lightly damp cloth or cotton swab.

10. Warranty

Implen guarantees that the product supplied has been thoroughly tested to ensure that it meets its published specification. The warranty included in the conditions of supply is valid for 24 months and is valid only if the product has been used according to the instructions supplied. Implen or your supplier can accept no liability of loss or damage arising from the faulty or incorrect use of this product.

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